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Instructions for Use on the Handling and Processing of Surgical Instruments

Our surgical instruments are designed and manufactured to meet the highest quality standards. We do not accept liability for products that have been modified, subjected to unintended use, or used improperly. Single-use products may not be reused.

The package inserts provided with the devices do not contain all of the information necessary for the selection and application of the devices. For proper handling, refer to other device-related instructions, such as the instructions on the surgical technique associated with the relevant system, as well as the special handling recommendations, where applicable.

Devices made of plastics (e.g. PA, PE, POM, UHMWPE) may not be located by means of an external imaging device.

1. Initial use

Non-sterile instruments must be cleaned, disinfected and sterilized in accordance with Section 4 prior to initial use. All users in central sterile supply departments and operating rooms should familiarize themselves with the function and compatibility of an instrument before it is used.

2. Performance test

Instruments must be tested for their fitness for operation and their functional efficiency every time before they are used. In addition, they must be devoid of visible residues and/or soiling. In particular, areas and functions such as measuring functions, compatibilities, shapes, cutting edges, tips, connections, closures, locks, detents, and moving components should be checked. Components and instruments made of plastics must be checked for wear due to ageing, e.g. cracks, embrittlement or chipping, and if such damage is evident they may no longer be used but must be replaced with new instruments.

The service life of instruments depends on materials, design, use, and processing. Therefore do not use damaged instruments or ones with evident changes in surfaces. Do not perform repairs yourself. Servicing and repairs may only be performed by the manufacturer's suitably qualified personnel.

3. Maintenance

After every use, subject the instruments to an inspection/performance test. After cleaning and thermal disinfection, hinge joints, screw threads and other moving parts must be coated with physiologically safe oil in order to reduce friction and wear. For instrument care we recommend: „neodisher IP Spray“ *). In individual cases also observe the additional, instrument-specific LINK cleaning and care instructions.

4. Cleaning and care

After every use, the instruments must be cleaned, disinfected and inspected by qualified personnel (e.g. DGSV competency). For this purpose they must be dismantled as far as possible. We recommend removing any heavy soiling from the instruments with a lint-free cloth immediately after use. The instruments must then be subjected to further cleaning before body fluids, tissue remnants and other residues are able to dry onto them. Therefore they should be sent to the central sterile supply department immediately after the operation.

4.1 Manual cleaning

- 4.1.1 Use the neutral-pH enzymatic cleaning solution prepared according to the manufacturer's instructions. The product we recommend is MediClean from Dr. WEIGERT, Hamburg.
- 4.1.2 Completely immerse the instrument (or its individual parts) in the enzymatic solution and allow it to soak for 20 minutes. Clean the instrument carefully using a brush with soft plastic bristles (taking special care with cavities and other poorly accessible areas) until all visible soiling has been removed. Lumens should be cleaned using a long, slender brush with soft bristles (e.g. a pipe cleaner). Note: The enzymatic solution should be changed if it has become heavily contaminated (with blood and/or turbidity).
Do not use metal brushes or steel wool for cleaning.
- 4.1.3 Take the instrument out of the enzymatic solution and for at least 3 minutes rinse it in purified water that has been obtained by ultrafiltration, RO (reverse osmosis), DI (deionization) and/or distillation, or a combination of them. Rinse out lumens, openings and other poorly accessible areas thoroughly.
- 4.1.4 Prepare the neutral-pH cleaning solution and pour it into an ultrasonic cleaner.
- 4.1.5 Completely immerse the instrument in the cleaning solution and clean ultrasonically for 10 minutes, preferably at 45–50 kHz.
- 4.1.6 Rinse the instrument thoroughly with purified water that has been obtained by ultrafiltration, RO, DI, and/or distillation, or a combination of them, either for at least 3 minutes or until the rinsing water contains no further signs of blood or soiling.
- 4.1.7 Repeat steps 4.1.5 and 4.1.6 using freshly prepared cleaning solution.
- 4.1.8 Dry the instrument with a clean, absorbent and lint-free disposable cloth.

4.2 Cleaning in a washer-disinfector

Before cleaning and disinfection in a validated washer-disinfector (WD) using qualified personnel (e.g. DGSV competency) we recommend:

- 4.2.1 Dry disposal in the operating room.
- 4.2.2 Complete disassembly of detachable components
- 4.2.3 Manual, mechanical removal of any heavy soiling, preferably under running water.
- 4.2.4 Ultrasonic cleaning in suitable washing baskets, e.g. in a BANDELIN RM 40/UH, with 400 mL SEPTO PRE CLEAN per 40 L municipal water at room temperature and 40 kHz for 15 minutes.

- 4.2.5 Lay the instruments in suitable washing baskets, either dismantled or opened. Care must be taken to ensure that the liquid can circulate round the instruments adequately and that there are no surfaces not being exposed to the liquid. Please observe the instructions provided by the manufacturer of the washer-disinfector. If instruments are hollow, this washing operation must be just as efficient internally as it is externally; use special-purpose washer-disinfector adapters where necessary.

- 4.2.6 Automatic cleaning and disinfection in a machine must be performed using a validated washer-disinfector, e.g. a BELIMED WD 290. For validated automatic cleaning (including thermal disinfection at 93 °C for 5 min) we recommend the following program parameters with an alkaline cleaning agent after steps 4.2.1 to 4.2.5 have been completed:

Program step	Water inlet	Step time	Metering time
Precleaning	Cold water	4 min	–
Cleaning	Demineralized water (DIW)	7 min	80 sec. at 40 °C
Neutralizing agent	DIW	1 min	20 sec.
Final rinse	DIW	5 min	10 sec. at 50 °C
Drying, 1st phase		5 min	–
Drying, 2nd phase		10 min	–

Program step	Temperature	Chemical
Precleaning	25 °C	–
Cleaning	60 °C	neodisher MediClean forte
Neutralizing agent	–	neodisherN
Final rinse	93 °C	neodisher Mediklar
Drying, 1st phase	120 °C	–
Drying, 2nd phase	120 °C	–

5. Sterilization

Our instruments are designed for sterilization in moist heat (steam sterilization) by the fractionated pre-vacuum method (134 °C and a holding time of at least 5 minutes), e.g. in a WEBECO EC 2120, operated by qualified personnel (e.g. DGSV competency).

Steam sterilization can be performed in standardized sterile paper packaging (e.g. DIN EN ISO 10993-12, 7405 and/or ANSI AAMI ISO 11607, 11737-2) or in a dedicated tray in a sterilization container conforming to e.g. DIN EN 868-8, DIN EN ISO 11607-1.

The sterilization method must be validated by the user. The temperature should not exceed 137 °C because otherwise handles, insulation or other non-metallic components might be damaged.

Care must be taken to ensure that before the instruments are sterilized they have been re-assembled and opened. With regard to handling and loading please follow the instructions provided by the manufacturer of your sterilizer.

Plastic components and plastic-coated components may not be sterilized in hot-air sterilizers. Recommendation: When using sterilization containers with trays, the filled trays should also be wrapped in suitable cloths or nonwovens. The total weight of each container should not exceed 10 kg.

If there is any suspicion of prion contamination, follow the recommendations issued by the Robert Koch Institute (www.rki.de) and keep to a holding time of at least 18 minutes at 134 °C for steam sterilization.

Further information about processing medical devices properly is available from the „Instrument Processing Working Group“ at www.a-k-i.org.

6. Handling and storage

Surgical instruments must always be treated with care, this particularly applies during transport, cleaning, maintenance, sterilization, and storage. The sterile status of the instruments depends, inter alia, on the sterile items packaging and the prevailing storage conditions and must be established together with the operator's hygiene officer on a case-by-case basis. Direct sunlight must be avoided. Improper handling and care, as well as unintended use, can lead to premature wear or damage.

7. Complaints about our products

Complaints of whatever kind must be filed with LINK Italia S.P.A. at:

complaint@linkitaliaspa.it. When filing a complaint, always quote the name or REF number of the relevant component along with the LOT number, your name and your contact address. Describe the reason for the complaint in brief.

8. Servicing

Medical devices that are sent in for servicing must be processed beforehand in such a way that they cannot constitute a hazard to third parties. We will be pleased to provide you with further information about special instrument sets, their applications, disassembly, cleaning, and care.

*) Products from Chemische Fabrik Dr. Weigert, Hamburg. For information about these products, amounts required and sources visit www.drweigert.com.

Please use only approved cleaning and disinfection chemicals listed on the DGHM list at www.vah-online.de.

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