-  **Patienteninformation zu Implantaten**
-  **Patient information for Implants**
-  **Информация за пациента за импланти**
-  **Informace o implantátech pro pacienty**
-  **Patientinformation angående implantater**
-  **Πληροφορίες προς τους ασθενείς για τα εμφυτεύματα**
-  **Información para el paciente sobre implantes**
-  **Patsienditeave implantaatide kohta**
-  **Implantteja koskeva potilastiedote**
-  **Informations aux patients relatives aux implants**
-  **Informacije o implantatima za pacijente**
-  **Betegtájékoztató implantátumokhoz**
-  **Informasi Pasien untuk Implan**
-  **Informazioni per il paziente sugli impianti**
-  **Informacija pacientui apie implantus**
-  **Pacienta informācija par implantiem**
-  **Patiëntinformatie voor implantaten**
-  **Pasientinformasjon om implantater**
-  **Informacja dla pacjentów o implantach**
-  **Informações para o paciente sobre implantes**
-  **Informații despre implanturi pentru pacienți**
-  **Информация об имплантатах для пациентов**
-  **Informácie pre pacientov v súvislosti s implantátmi**
-  **Informacije za paciente za vsadke**
-  **Information om implantat till patienter**
-  **İmplantlar için Hasta Bilgisi**
-  **Інформація для пацієнтів стосовно імплантатів**
-  **植入体患者的注意事项**

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1. Dear patient,

This document contains important information about your joint replacement. It shall help you to get fit again and to stay mobile.

2. Implant Lifetime

The lifetime of our implants is limited in principle and is determined by individual factors such as, for example, body weight and the level of activity of the patient, as well as by the quality and professional execution of the implantation. Based on these individual influencing factors, Waldemar Link defines the overall average lifetime of an implant based on its survival rate (i.e. the proportion of functional implants after a certain period of time starting from the time of implantation). According to the results of the tests performed, the survival rate of our implants corresponds to the general state of the art at the time of approval of the implants.

3. CMR

Some system components contain cobalt as an alloy ingredient in a concentration above 0.1 % weight by weight. Cobalt is listed as a hazard substance (CMR) and assigned to the Hazard Class and Category Code: Carc. 1B, Repr. 1B.

4. Important Information

All joint replacement surgeries

- Go to all follow-up examinations and contact the doctor, if problems occur.
- An implant card is given to all patients who have received an implant from a healthcare facility or clinic. The implant card contains all the necessary information about the implant. In addition, it offers those affected the opportunity to identify themselves at airports or at other security checks with metal detectors, thus avoiding problems.

Tips for everyday life after a hip replacement surgery

- Avoid heavy physical work.
- Eliminate tripping hazards in your living area, such as loose carpet edges.
- Avoid sports with heavy impact or a high risk of injury.
- Suitable sports are light gymnastics, hiking on flat terrain and cycling in areas with only slight inclines.
- Swimming is also allowed with a hip endoprosthesis, preferably in crawl style.

Tips for everyday life after a knee replacement surgery

- Eliminate tripping hazards in your living area, such as loose carpet edges.
- Use a backpack to shop, making heavy loads easier.
- Avoid sports with heavy impact and abrupt changes of direction such as jogging, soccer, tennis and horse riding.
- Suitable sports are cycling, gymnastics, hiking and swimming.
- It is best to choose sports that you are already familiar with.

Tips for everyday life after a shoulder replacement surgery

- Avoid heavy physical work.
- Avoid sports with heavy impact or a high risk of injury.
- Prefer sports with low shoulder and arm strain and low risk of falling.
- Suitable sports are aqua jogging, aerobics, gymnastics, jogging, canoeing, moderate weight training under supervision or instruction, swimming, cycling, cross-country skiing, dancing, hiking and Nordic walking.
- It is best to choose sports that you are already familiar with.
- Consult with your treating physician regarding the types of sports that are appropriate for you.

Tips for everyday life after a revision surgery

- Avoid heavy physical work.
- Eliminate tripping hazards in your living area, such as loose carpet edges.
- Keep a healthy body weight.
- Exercise regularly adapted sports. Swimming is very suitable, preferably in crawl style.
- Do not engage in sports that involve strong impacts or a high risk of injury, such as tennis or squash.
- Consult with your treating physician, if you have questions about any sport.

Tips for your next flight or other security checks

- You have a joint replacement implanted, partly consisting of metal. This may interfere with metal detection devices. Be prepared to show your implant card. The implant card helps the staff to determine the cause of detector activation.

Explanation of Pictograms

	Manufacturer
	Patient information website
	Material (number)
	CE-Symbol (Certificated product)
	Contains hazardous substances

Legend of Materials

	
M-No.	Material description
M1	Cobalt-based alloy, CoCrMo, ISO 5832-12; ASTM F1537, (EndoDur-S) Co 63-69% Cr 26-30% Mo 5-7% Ni, Fe, C, Mn, Si, N < 1%
M2	Titanium-aluminum alloy, Ti6Al4V, DIN EN ISO 5832-3; ASTM F136, (Tilastan-S) Ti 89-91% Al 5.5-6.5% V 3.5-4.5% Fe, C, H, N, O < 0.25%
M3	Ultra-high molecular weight polyethylene, UHMWPE, ISO 5834-2 / ASTM F-648 ASTM F-648 UHMWPE 100%
M4	Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene, UHMWPE, ISO 5834-2 / ASTM F-648 / ASTM F-2565, (X-Linked PE) UHMWPE 100%

M5	Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene, UHMWPE, with vitamin E, ISO 5834-2 / ASTM F-648 / ASTM F-2565 / ASTM F-2695, (E-DUR) UHMWPE 99.9% α -Tocopherol 0.1%
M6	Calcium phosphate coating, ASTM F-1609, (HX) CaP 100%
M7	Titanium nitride coating, TiNbN, titanium niobium nitride, (LINK PorEx) N 50 At% $\pm$ 20% Ti 42 At% $\pm$ 20% Nb 8 At% $\pm$ 20% O, C, Fe < 5At%
M8	Polyphenylsulfone, PPSU
M9	Neodymium, NdFeB, Nickel-plated
M10	Stainless steel, X5CrNi18-10, DIN EN 10088-1
M11	Stainless steel, X8CrNiS 18-9, DIN EN 10088-1
M12	Stainless steel, X20Cr13, DIN EN 10088-1
M13	Stainless steel, X10CrNi 18-8, DIN EN 10270-3
M14	Stainless steel, X90CrMoV18, DIN EN 10088-3
M15	Stainless steel, X5CrNiCuNb 16-4 / X5CrNiCuNb 17-4, ISO 7153-1
M16	Stainless steel, X2CrNiTiMo12-11-2-1
M17	Stainless steel, X17CrNi16-2, DIN EN 10088-3
M18	Stainless steel, XM-16, ASTM F 899-12
M19	Cobalt-based alloy, CoCrNiMoFe, ISO 5832-7 ASTM F1058
M20	Aluminium alloy, Al99,5, EN AW-1050A, DIN EN 573-3
M21	Silicone, BGA XV
M22	Polyamide, PA6 SA, DIN EN 15860
M23	Polypropylene homopolymer, PP-H, DIN EN 15860, ASTM D4101
M24	Polyvinyl chloride, PVC
M25	Polytetrafluoroethylene. PTFE
M26	Stainless steel, X5CrNiMo17-12-2, DIN EN 10088-3, AISI 316
M27	Stainless steel, X6CrNiMoTi17-12-2, DIN EN 10088-2
M28	Stainless steel, X46 Cr-13, ISO 7153-1
M29	Stainless steel, X2CrNiMo 18-15-3, DIN ISO 5832-1:97, ASTM F138-00, ASTM F139-00
M30	Cobalt-based alloy, CoCrMo, DIN ISO 5832-4, ASTM F75 (EndoDur)
M31	Titanium-aluminum alloy, Ti6Al4V, ASTM F1108, (Tilastan)
M32	Titanium-aluminum alloy, Ti6Al4V, DIN EN ISO 5832-3, (Tilastan-E)
M33	Stainless steel, X2CrNiMo 17-12-2, ISO 7153-1
M34	Stainless steel, X2CrNiMo18-14-3, DIN EN 10088-3
M35	Aluminium alloy, AlMg1, EN EN5005/H14, DIN EN 573-3
M36	Thermoplastic elastomer, TPE
M37	Silicone rubber, FDA§ 177.2600
M38	Cobalt-based alloy, CoNi35Cr20Mo10, ISO 5832-6
M39	Titanium-aluminum alloy, Ti6Al4V, ASTM F2924
M40	Commercially Pure Titanium coating, CPTi, ASTM F1580
M41	Aluminium alloy, AlSi1MgMn, DIN EN 573-3
M42	Stainless steel, X105CrMo17, DIN EN 10088-3; EN 10278
M43	Zirconia toughened aluminium, ZTA, ISO 6474-2
M44	Polyetheretherketone, XL 515 CF, PEEK

M45	Silicone Rubber
M46	Stainless steel, FeCrNiMnMoNbN, ISO 5832-9, ASTM F1586-13
M47	Titanium-aluminum alloy, Ti6Al4V, DIN EN ISO 5832-3
M48	Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene, UHMWPE, ISO 5834-1/2 / ASTM F-648
M49	Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene with vitamin E, UHMWPE, ISO 5834-1/2 / ASTM F-648 / ASTM F-2695
M50	Cobalt-based alloy, CoCrMo, ISO 5832-12; ASTM F799
M51	Stainless steel, X20CrNiMoS13-1, ISO 7153-1
M52	Polysulfone, PSU
M53	Stainless steel, X40CrMoVN16.2, ASTM F 899
M54	Ethylene propylene diene monomer rubber, EPDM rubber
M55	Polyamide, PA12, DIN EN 15860
M56	Cast steel, GX5CrNiMoNb19-11-2, DIN EN 10283mod
M57	Cast steel, GX5CrNiMoNb19-11-2, DIN EN 10283
M58	Stainless steel, X2CrNi18-9, DIN EN 10088-1
M59	Silicone Rubber, RTV-4136-M
M60	Stainless steel, X2CrNiMo 18-15-3, DIN ISO 5832-1:97, ASTM F138-00, ASTM F139-00
M61	Stainless steel, X38CrMoV15, DIN EN ISO7153-1
M62	Stainless steel, X30Cr13, DIN EN ISO7153-1
M63	Polyoxymethylene, POM-C
M64	Stainless steel, S21800, ASTM F899-12b
M65	Commercially Pure Titanium coating, CpTi, DIN ISO 5832-2, ASTM F1580, (TiCaP)
M66	Stainless steel, X40 CoorVN 16-2, DIN EN 10088-3, AISI 420 MOD, EN 10278-h8
M67	Silicone Rubber, VMQ
M68	Polyamide, PA 2200
M69	Stainless steel, X15Cr13, AISI 410D, ISO 7153-1
M70	Aluminum Oxide, Al ₂ O ₃ , ISO 6474-1
M71	Polyethylene, HDPE
M72	Polymethyl methacrylate, PMMA
M79	Stainless steel, X40 CrMoVN 16-2, DIN EN 10088-3, AISI 420 MOD, EN 10278-h8
M80	Stainless steel, X4CrNi 18 12, ISO 7153-1
M81	Stainless steel, X10Cr15, AISI 429
M82	Stainless steel, X90CrMoV17, DIN EN 10088-3
M83	Aluminium alloy, AlMg3
M84	Stainless steel, X3CrNiCu18-9-4, AISI 304 Cu
M85	Stainless steel, X39CrMo17-1, DIN EN 10088-3; EN 17440; ISO7153-1

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, HX, TiCaP, X-LINKed, PorAg, LINK PorEx, BiPorEx, PorEx-Z, TrabecuLink, Tilastan, customLink, RescueSleeve, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

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